

# TopCount

TCA-028

## Dual Label DPM on the TopCount™

### Introduction

In dual label assays, the two isotopes are detected at different efficiencies due to the different energies of their beta emissions. Because the energy of a beta particle is distributed from 0 to  $E_{\max}$ , overlap of the two isotope spectra always occurs. DPM correction algorithms subtract the contribution caused by this overlap (spill) and adjust for the relative efficiencies. The use of DPM allows for meaningful interpretation and comparison of results in a dual label assay.

There are two methods for performing dual label DPM on the TopCount. The first is for samples that have variable levels of quench. This method requires the use of an external standard and the creation of two single label quench standard sets for DPM correction. The external standard option is available with the 24-well plate format in the TopCount. The second method of dual label counting assumes that all the samples are prepared in the same manner, and have a constant level of quench. Efficiency values

from known standards are used to correct CPM to DPM. This constant quench method can be used for samples in either the 24- or the 96-well format. This paper will discuss the constant quench method in detail. Variable quench dual label counting performance is presented in TopCount Topics TCA-015, "Quench and Quench Correction."

### Constant Quench Dual Label DPM

The first step in setting up constant quench dual label counting is to determine the two appropriate counting regions. The lower energy isotope should be fully detected in the lower energy region. The region for the higher energy isotope should be set to exclude any counts from the lower energy isotope; that is, to exclude any spill up (see Figure 1). RegionView is a useful feature on the TopCount for determining the appropriate region settings. (Refer to the TopCount Reference Manual.) Once the settings are known, a dual label nuclide should be created in the Nuclide Library.




Figure 1.

Example: For  $^3\text{H}/^{14}\text{C}$  using MicroScint™ cocktail, version 4.0 or higher TopCount software.

### Create a Custom Nuclide Setting

(See Figure 2)

From the STATUS window, press F1 - Directory, then highlight protocol of interest and press F3 - Count Conditions

Highlight the Radionuclide? line, and press F3 - Library

Select the function key F3 - Add Nuc to add a new nuclide setting

Type ' $^3\text{H}/^{14}\text{C}$  Liquid' <enter>

Select 'Liq/Plast' scintillator, 'Low' energy, and 'Normal' counting mode

Region A: LL: type '2.9' <enter>, UL: '30' <enter>

Region B: LL: type '35' <enter>, UL: '256' <enter>

Press F1 to Exit/Save new settings

After defining the counting regions, individual standards of each isotope must be counted under these regions to determine the efficiencies, and spill factors for use in correcting CPM to DPM. The standards should be of known DPM and prepared under the same final chemical and physical conditions (same quench level) as the assay unknowns. After counting the standards, calculate the efficiency for the low energy isotope in Region A, the efficiency of

the high energy isotope in Region B, and the spill. Spill is defined as the amount of high energy isotope that is measured in Region A.

The equations for these calculations are as follows:

$$1. \text{ } ^3\text{H Efficiency} = \frac{\text{CPM}_A}{\text{DPM}_{3\text{H}}} = \text{EFF}_{\text{LA}}$$

$$2. \text{ Spill Factor} = \frac{\text{CPM}_A}{\text{DPM}_{14\text{C}}} = \text{EFF}_{\text{HA}}$$

$$3. \text{ } ^{14}\text{C Efficiency} = \frac{\text{CPM}_B}{\text{DPM}_{14\text{C}}} = \text{EFF}_{\text{HB}}$$

Note: LA = Low energy isotope in Region A  
 HB = High energy isotope in Region B  
 EFF<sub>HA</sub> = Spill factor for high energy isotope into Region A

Using the values from the above equations, custom cells in the TopCount can be created to automatically provide the corrected DPM values on the disk file and printout. The equations for DPM (broken down into steps) are as follows:

1.  $\text{DPM}_{14\text{C}} = \text{CPM}_B / \text{EFF}_{\text{HB}}$
2.  $\text{SPILL} = \text{DPM}_{14\text{C}} * \text{EFF}_{\text{HA}}$
3.  $\text{CCPM}^\dagger_{3\text{H}} = \text{CPM}_A - \text{SPILL}$
4.  $\text{DPM}_{3\text{H}} = \text{CCPM}_{3\text{H}} / \text{EFF}_{\text{LA}}$

† CCPM = Corrected CPM

aF1-Help aF2-Count Ctrl aF3-Decay aF4-RegionView aF5-Move Window aF6-Mode						
Protoc			Prot# Protocol Name			
3H/14C Liquid Nuclide						
Sample	▲	Energy	Efficiency	Regions		
Nuclide	Scintillator	Range	Mode	A:LL - UL	B:LL - UL	
3H/14C Liquid	Liq/Plast	Low	Normal	2.90 - 30.00	35.00 - 256.0	
Enter application-specific nuclide name.						
F1-Exit/Save		F2-Cancel		F3-Clear		
125I-VS-SPA	Glass	Low	Hi Sen	0.00 - 100.0	0.00 - 256.0	
3H/14C-MicrSnt	Liq/Plast	Low	Normal	2.90 - 10.00	10.00 - 256.0	
32P-Cerenkov	Liq/Plast	Low	Hi Eff	7.50 - 256.0	7.50 - 256.0	
32P-MicroScint	Liq/Plast	High	Normal	2.90 - 150.0	2.90 - 150.0	
45Ca-MicroScnt	Liq/Plast	High	Normal	2.90 - 150.0	2.90 - 150.0	
NEW_45Ca-Micro	Liq/Plast	High	Normal	2.90 - 150.0	2.90 - 150.0	
NEW_51Cr-Flash	Liq/Plast	Low	Normal	2.90 - 150.0	2.90 - 256.0	
NEW_3H-MicroSc	Liq/Plast	Low	Normal	2.90 - 35.00	2.90 - 256.0	
3H/14C Liquid	Liq/Plast	Low	Normal	2.90 - 30.00	35.00 - 256.0	
Use Up or Down arrow keys to select an entry						

Figure 2.

These equations are then entered into the custom cells for the protocol of interest.

### Custom Cell Definition

(See Figure 3)

From the Directory page, highlight the protocol of interest and press F8 - Format Data Cells. This opens the data cell edit window. (Use arrow keys to move the cursor in this window.)

- Scroll down past the preset cells and highlight the first custom cell (C# 20)  
Type "BXXXXXX" <enter> for format  
Highlight the name 'Custom 1' and type "DPM-14C" <enter>
- Type "C7/EFF<sub>HB</sub>" <enter>; where EFF<sub>HB</sub> is the decimal value determined above
- Move the cursor to the next custom cell  
Type "BXXXXXX" <enter> for format  
Highlight the name 'Custom 2' and type "SPILL" <enter>

Type "C20\*EFF<sub>HA</sub>" <enter>; where EFF<sub>HA</sub> is the decimal value determined above

- Move the cursor to the next custom cell  
Type "BXXXXXX" <enter> for format  
Highlight the name 'Custom 3' and type "C-CPMA" <enter>  
Type "C6-C21" <enter>
- Move the cursor to the next custom cell  
Type "BXXXXXX" <enter> for format  
Highlight the name 'Custom 4' and type "DPM-3H" <enter>  
Type "C22/EFF<sub>LA</sub>" <enter>; where EFF<sub>LA</sub> is the decimal value determined above  
Press F1 - Exit to save changes and return to the Directory page

The cells containing the final result for DPM<sub>3H</sub> and DPM<sub>14C</sub> must be selected in the 'Printer Output' window (Figure 4) and the 'Disk File Output' window (Figure 5) to get the calculated values in the file or on the hard copy respectively.

C#	Name	Format	Name	Equation
16	Type	BXXXX		
17	FLAG	BXXXX		
18	{spaces}	BBBBBBBB		
19	{CR-LF}	X		
20	DPM-14C	BXXXXXX	DPM-14C	C7/0.6
21	SPILL	BXXXXXX	SPILL	C20*0.15
22	C-CPMA	BXXXXXX	C-CPMA	C6-C21
23	DPM-3H	BXXXXXX	DPM-3H	C22/0.33
24	Custom 5	BXXX.XXX	Custom 5	
25	Custom 6	BXXX.XXX	Custom 6	

Figure 3.

**Example: Constant Quench Dual Label DPM**

The single label standards are counted using the nuclide setting described previously. Each standard is known to contain 100,000 DPM.

For the <sup>3</sup>H standard:

$$CPM_A = 33,000 \text{ and } CPM_B = \text{Background}$$

For the <sup>14</sup>C standard:

$$CPM_A = 15,000 \text{ and } CPM_B = 60,000$$

The efficiency and spill values are then calculated:

$$E_{HB} = 60,000/100,000 = 0.6 \text{ (60\% efficiency)}$$

$$E_{LA} = 33,000/100,000 = 0.33 \text{ (33\% efficiency)}$$

$$E_{HA} = 15,000/100,000 = 0.15 \text{ (15\% efficiency)}$$

An unknown sample containing <sup>3</sup>H and <sup>14</sup>C labeled material is counted and the results are as follows:

$$CPM_A = 10,000 \text{ } CPM_B = 10,000$$

Are there equal amounts of both <sup>3</sup>H and <sup>14</sup>C in the sample?

$$DPM_{14C} = 10,000/0.6 = 16,667 \text{ DPM of } ^{14}C$$

$$SPILL = 16,667 * 0.15 = 2,500$$

$$CCPM_{3H} = 10,000 - 2,500 = 7,500 \text{ CPM of } ^3H$$

$$DPM_{3H} = 7,500/0.33 = 22,727 \text{ DPM of } ^3H$$

**Conclusion**

In reality, there is almost 1.4 times as much <sup>3</sup>H present as <sup>14</sup>C in this sample. Using DPM allows for an accurate comparison of the amounts of the two isotopes in a given sample, and for comparison of values from different experiments. The custom cell feature on the TopCount provides a convenient, automatic means of performing constant quench dual label assays with high throughput for screening applications.




Figure 4.




Figure 5.