SNP Genotyping with Fluorescence Polarization Detection

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
Figure 1. Fluorescence polarization (FP). The FP of a molecule is proportional to the molecule's rotational relaxation time, which is a function of the solvent viscosity, absolute temperature, and molecular volume. If the viscosity and temperature are held constant, FP is directly proportional to the molecular volume, which is directly proportional to the molecular weight. If the fluorescent molecule is large, it tumbles more slowly in space and FP is preserved. If the molecule is small, it tumbles faster and depolarized emission is observed.

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
Figure 3. FP analysis of primer extension assay results. The homozygous C samples have high FP values for TAMRA but low FP values for R6G and occupy the lower right corner of the plot. The homozygous T samples have high FP values for R6G but low FP values for TAMRA and occupy the upper left corner of the plot. The heterozygotes are found in the upper right corner and have high FP values for both dyes. The negative controls are found near the origin, with low FP values for both dyes.

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
Figure 4. 5'-Nuclease assay with FP detection.

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
Figure 5. Fluorescence polarization analysis of 5'-nuclease assay with marker D2S1301 on 90 individuals. The homozygous allele 1 samples have low VIC FP values but high FAM FP values and occupy the lower right corner. The homozygous allele 2 samples have low FAM FP values but high VIC FP values and occupy the upper left corner. The heterozygotes have low FP values for both dyes and are found near the origin. The negative controls occupy the upper right corner, with high FP values for both dyes.

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
Figure 6. Invader Assay. When there is a match between allele-specific flap probe (primary probe) and DNA target at the polymorphic site, an overlapping structure between the 1° Invader probe and the primary probe is formed. A thermostable Cleavase enzyme recognizes this structure and cleaves the allele-specific probe, which releases the 5' flap. This 5' flap in turn serves as the 2° Invader probe in the second reaction where the signal probe was cleaved and the fluorescent molecule is released.

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
Figure 7. PCR-Invader genotyping data for TSC 1469. Homozygous T individuals form a cluster in the upper left corner with high Tamra FP values and low Fluorescein FP values. Homozygous G samples form a cluster in the lower right corner with high Fluorescein FP values and low Tamra FP values. Heterozygotes have low Tamra and Fluorescein FP values, and they cluster in the lower left corner. The no-DNA negative controls, with both signal probes intact, have high FP values for both dyes analyzed and occupy the right upper corner of the plot, indicating that both dyes are part of high molecular weight species.

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References

1. Afonina I, Zivarts M, Kutyavin I, Lukhtanov E, Gamper H, Meyer RB. 1997. Efficient p priming of PCR with short oligonucleotides conjugated to a minor groove binder. *Nucleic Acids Res* 25:2657-2660.
2. Checovich WJ, Bolger RE, Burke T. 1995. Fluorescence polarization a new tool for cell and molecular biology. *Nature* 375:254-256.
3. Chen X, Levine L, Kwok PY. 1999. Fluorescence polarization in homogeneous nucleic acid analysis. *Genome Res* 9:492-498.
4. Gibson NJ, Gillard HL, Whitcombe D, Ferrie RM, Newton CR, Little S. 1997. A homogeneous method for genotyping with fluorescence polarization. *Clin Chem* 43: 1336-1341.
5. Heyduk T, Ma Y, Tang H, Ebright RH. 1996. Fluorescence anisotropy: rapid, quantitative assay for protein-DNA and protein-protein interaction. *Methods Enzymol* 274:492-503.
6. Hsu TM, Law SM, Duan S, Neri BP, Kwok PY. 2001a. Genotyping single nucleotide polymorphisms by the Invader assay with dual-color fluorescence polarization detection. *Clin Chem* 47:1373-1377.
7. Hsu TM, Chen X, Duan S, Miller RD, Kwok PY. 2001b. A universal SNP genotyping assay with fluorescence polarization detection. *Biotechniques* 31:560-570.
8. Kaiser MW, Lyamicheva N, Ma W, Miller C, Neri B, Fors L, Lyamichev VI. 1999. A comparison of eubacterial and archaeal structure-specific 5'-exonucleases. *J Biol Chem* 274:21387-21394.
9. Kutyavin IV, Lukhtanov EA, Gamper HB, Meyer RB. 1997. Oligonucleotides with conjugated dihydropyrroloindole tripeptides: base composition and backbone effects on

hybridization. *Nucleic Acids Res* 25:3718-3723.


10. Kwok PY. 2001. Methods for genotyping single nucleotide polymorphisms. *Annu Rev Genomics Hum Genet* 2:235-258.

11. Latif S, Bauer-Sardina I, Ranade K, Livak KJ, Kwok PY. 2001. Fluorescence polarization in homogeneous nucleic acid analysis II: 5'-nuclease assay. *Genome Res* 11: 436-440.

12. Lee LG, Livak KJ, Mullah B, Graham RJ, Vinayak RS, Woudenberg TM. 1999. Seven-color, homogeneous detection of six PCR products. *Biotechniques* 27:342-349.

13. Livak KJ. 1999. Allelic discrimination using fluorogenic probes and the 5' nuclease assay. *Genet Anal* 14:143-149.

NEED A MICROPLATE READER?

- RAPID KINETIC ASSAYS
- CELL BASED ASSAYS
- ENZYME ACTIVITY
- REPORTER GENE
- DNA/RNA QUANTIFICATION
- SECONDARY MESSENGERS
- FRET
- IMMUNOASSAYS
- ROBOTIC COMPATIBLE
- ELISA

FLUORESCENCE POLARIZATION
 LUMINESCENCE INTENSITY

 BMG Labtechnologies
www.bmglabtech.com